


Get better with ActivaPin™

Easy insertion with pin applicator


Since 2007

BETTER FIXATION

- For fixation of bone fractures, osteotomies, arthrodesis and osteochondral fractures in the presence of appropriate immobilization.
- ActivaPin™ is an implant with ultra-high-strength properties and it has an excellent track record (according the US Maude database). More than 80,000 implants have been used worldwide and the occurrence of complications is less than 0.01%.
- ActivaPin™ decreases stress shielding in bones and biodegrades in the body safely and in a controlled manner.



An instant locking effect into the drill hole. Patented grooved surface design to provide enhanced rotational stability.



Example of an indication:
Chevron osteotomy
with ActivaPin™

**FDA cleared and
CE approved**

SAFE TO USE


- Bending modulus is close to the value of cortical bone, contrary to metals.
- ActivaPin™ maintains its intended function for at least 8 weeks.
- Bioabsorption reduces the risk of long-term complications.
- Aseptic surgical procedure due to packaging, instruments and implants.
- Implants have an excellent clinical track record relating to surgical use worldwide.

EXCELLENT PRODUCT PROPERTIES

- High-strength properties offer easy insertion and safe medical use.
- Complete bioabsorption takes place within approximately 2 years, eliminating the need for implant removal surgery.
- Implant is supplied gamma sterilized.
- Designed and manufactured in Finland according to certified quality management system.

J.B. Implants Services bv

Robijnborch 7 T 073 52 20780
5241 LK Rosmalen F 073 52 21245
www.jb-implants.nl E info@jb-implants.nl



bioretec

Bioretec develops and produces globally innovative bioabsorbable products and surgical solutions for better patient outcome, user friendliness and cost-efficiency in clinical care.

Bioretec Ltd. • Hermiankatu 22 • FI-33720 Tampere • Finland
Tel. +358 20 778 9500 • Trade Reg. No. 738.932 • VAT No. FI14741969
sales@bioretec.com • www.bioretec.com

