

SCRU™ 2©

True headless screw

- Cannulated over 1mm KW
- Compressive
- Dual pitch
- Self drilling
- Self tapping



Case from: Michel Merle, M.D.

AREX®
Stay the course



SCRU™ 2© LOW PROFILE

2.5 mm $\varnothing$ top of the first thread
2.3 mm $\varnothing$ top of the second thread



References


Length	Low Profile
20 mm	P2P20 LP
25 mm	P2P25 LP
30 mm	P2P30 LP
35 mm	P2P35 LP

SCREW PENETRATING TIP
Self Drilling - Self Taping

SECOND THREAD
Self Compressive

LP: 2.5 mm $\varnothing$ top of the thread

LP: 2.3 mm $\varnothing$ top of the thread



CE 0499

Manufacturer:

AREX®

BP20

5 allée du clos tonnerre

F91125 Palaiseau cedex FRANCE

www.arex.fr

J.B. Implants Services bv

Robijnborch 7
5241 LK Rosmalen
www.jb-implants.nl

T 073 52 20780
F 073 52 21245
E info@jb-implants.nl